
Innovative Products for Life Sciences

AquaClean

Microbiocidal Additive for Water Soluble Heating Bath Fluids





AquaClean inhibits growth of bacteria, algae, fungi and precipitation of inorganic salts.

AquaClean shows a blue indicator color. Bacterial contamination is indicated by decoloration of the fluid.

AquaClean does not contain aldehydes or azides.

Application:

- Add 5ml AquaClean to 1 liter of heating bath fluid.
- Renew fluid and additive every 4 weeks.
 In case the blue color of fluid has faded, renew earlier.

Compatible for use as an additive in water reservoirs inside cell culture incubators. (Its non-volatile active ingredients cannot invade cell-cultures via air.)

PL-WAK-aquaclean-eng-01-r01-05/08

Cat.#	WAK-AQ-250-50L 250 ml concentrate sufficient for 50 liters of fluid	Active ingredients:	Surfactants,quaternary ammonium compounds, complexing agents
Precautions:	Avoid eyes and skin contact	Storage:	Room temperature
Disposal:	Biologically degradable	Materials compatibility:	Universal except for aluminium and acrylic glass

Siemensstraße 9 D-61449 Steinbach/Ts. • Germany Tel. +49 61 71 / 28 43 0 • Fax +49 61 71 / 69 56 46 e-mail: info@wak-chemie.com • www.wak-chemie.com

WAK - Chemie Medical GmbH



		Material Safety Data Sheet	
1.	Product Description Catalog No.: Product name: Purpose of use:	WAK-AQ-250-50L AquaClean Microbiocidal additive for water soluble heating bath fluids.	
2.	Chemical ingredients Complexing agents, quaternary ammonium compounds, surfactants, demineralized water.		
3.	Eye or mucosa contact: rin	sh skin with water and mild soap. se thoroughly with water for 15 minutes. Contact a physician. ve clothing to avoid body contact.	
4.	Emergency and first aid procedures In case of ingestion: rinse mouth with water. Contact a physician. Never give anything by mouth to an unconscious person. Never make an unconscious person vomit. In case of allergic reactions contact a physician.		
5.	Fire and explosion hazard data Non inflammable. Water-soluble.		
6.	Steps to be taken in case material is released or spilled Use absorbent material for removal. Thereafter, rinse with water.		
7.	Handling instructions Avoid contact with eyes or mucosa. Wear protective clothing, gloves and glasses.		
8.	Storage Store light protected at room temperature.		
9.	Physical data Boiling point: Appearance: pH-value: Odor: Specific weight: Solubility: Flammable limits: Chemical stability:	100°C blue liquid 7.1 n/a 1.01 g/cm ³ water soluble none stable	
10.	Stability data Stable at ambient temperature. No decomposition by boiling or evaporating.		
11.	Toxicological data Potentially toxic if ingested.		
12.	Waste disposal methods Diluted 1:200 liquid can be disposed like grey water .		
13.	Transportation requirements No special instructions.		
14.	General While the information and recommendations set forth herein are believed to be accurate as of the date hereof, WAK-Chemie Medical GmbH makes no warranty with respect thereto and disclaims all liability from reliance thereon. Last revised 28. March 1999		

03/99-ConTaFree-1017-AC-SDB-E

Siemensstraße 9 D-61449 Steinbach / Ts. Germany Tel. + 49 61 71 / 2 84 30 • Fax + 49 61 71 / 69 56 46 e-mail: info@wak-chemie.com • www.wak-chemie.com

WAK - Chemie Medical GmbH

