# **Innovative Products for Life Sciences**

# **Biocidal ZF**<sup>™</sup>

Spray Disinfectant for Incubators and Sterile Cabinets in Cell Culture Area





# Effectiveness:

Biocidal ZF<sup>™</sup> is effective against bacteria, spores and enveloped viruses.

## Non-volatile:

Active microbiocidal ingredients are non-volatile. They protect cell cultures from microbial contamination and do not invade cell cultures via air. Thus cell cultures are protected against contamination and the disinfectant itself.

## **Application for incubators**

- Empty incubator
- Spray Biocidal ZF<sup>™</sup> on all surfaces and let settle
- Put cultures right back

Surfaces will be protected for 2-4 weeks. (Before disinfecting you may clean incubator with Biocidal ZF™ if necessary)

### **Application in sterile benches**

- Empty sterile bench
- Then clean with Biocidal ZF<sup>™</sup> and finally spray on all surfaces.

Protection will prevail until next workshift. No combustible gases can form, unlike by use of ethanol.

Cat.#		1 x 1 6 x 1 12 x 0.2	liter spray bottle liter spray bottle liter spray bottle	Active Ingredients: Storage:	Quaternary ammonium compounds -free of aldehydes- Room temperature
Precautions: Avoid eyes and skin contact				Materials compatibility:	Universal inclusive fittings
Disposal:	Biologically degradable				within incubators and sterile cabinets

Siemensstraße 9 D-61449 Steinbach / Ts. Germany Tel. + 49 61 71 / 2 84 30 • Fax + 49 61 71 / 69 56 46 e-mail: info@wak-chemie.com • www.wak-chemie.com

# WAK - Chemie Medical GmbH



#### **Material Safety Data Sheet**

#### 1. Product Description

Catalog No.: Product name: Purpose of use: WAK-ZF-1; WAK-ZF-6; WAK-ZF-12 Biocidal ZF™ Ready-to-use disinfectant to protect cell cultures from microbial contamination.

#### 2. Chemical ingredients

Quaternary ammonium compounds, auxiliary additives.

#### 3. Health hazard data

Weakly irritant for skin. Irritant for eyes and mucosa. Do not inhale spray-aerosols. Skin contact: promptly wash skin with water and mild soap.Eye or mucosa contact: rinse thoroughly with water for 15 minutes. Consult physician. Wear appropriate protective clothing to avoid body contact.

#### 4. Emergency and first aid procedures

Ingestion: Give large amounts of water. Never give anything by mouth to an unconscious person. Never try to make an unconscious person vomit. Contact a physician.

#### 5. Fire and explosion hazard data

Non inflammable. Water-soluble.

#### 6. Steps to be taken in case material is released or spilled

Take up with absorbent material, rinse with water.

#### 7. Handling instructions

Avoid contacts with eyes or mucosa. Wear protective clothing, gloves and glasses.

#### 8. Storage

Protect from mechanical damage. Store light protected at room temperature

#### 9. Physical data

Boiling point:	100°C
Appearance:	colorless liquid
pH-value:	8.5 ± 0.5
Odor:	n/a
Specific weight:	$0.99 \pm 0.01 \text{ g/cm}^3$
Solubility:	water soluble
Flammable limits:	none
Chemical stability:	stable > pH6

#### 10. Stability data

Stable. Dangerous decomposition products not known.

#### 11. Toxicological data

Potentially hazardous if inhaled or ingested.  $LD_{so}$  oral rat 50 000 mg/kg.

#### 12. Waste disposal methods

Biologically degradable. Dispose like grey water.

## **13.** Transportation requirements

No special instructions

#### 14. General

While the information and recommendations set forth herein are believed to be accurate as of the date hereof, WAK-Chemie Medical GmbH makes no warranty with respect thereto and disclaims all liability from reliance thereon. Last revised 12- September 2007

#### 09/07-ConTaFree-1016B-SDB-E

Siemensstraße 9 D-61449 Steinbach/Ts. • Germany Tel. +49 61 71 / 28 43 0 • Fax +49 61 71 / 69 56 46 e-mail: info@wak-chemie.com • www.wak-chemie.com

# WAK - Chemie Medical GmbH

