#### **Innovative Products for Life Sciences**

# Acryl AquaClean

Microbiological Additive for Heating Bath Fluids Compatible with Acrylic Glass and Aluminium

## Water Baths



biologically degrad

Acryl AquaClean inhibits growth of bacteria, algae and fungi.

Acryl AquaClean is long-term compatible with equipment containing acrylic glass and aluminium materials.

Acryl AquaClean does not contain azides or aldehydes.

#### **Application:**

Add 5ml AquaClean to 1 liter of heating bath fluid.

Renew fluid and additive every 4 weeks. In case the blue color of fluid has faded, renew earlier.

Compatible for use as an additive in water reservoirs inside cell culture incubators. (Its non-volatile active ingredients cannot invade cell-cultures via air.)

| Cat.#        | WAK-AQA-250-50L<br>250 ml concentrate<br>sufficient for 50 liters of fluid | Active Ingredients:<br>Storage: | Ketones, Surfactants<br>Room temperature |
|--------------|--|---------------------------------|--|
| Precautions: | Avoid eyes and skin contact  | Materials compatibili           | ity: Universal                           |
| Disposal:    | Biologically degradable  |                                 |  |

Siemensstraße 9 D-61449 Steinbach / Ts. Germany Tel. + 49 61 71 / 2 84 30 • Fax + 49 61 71 / 69 56 46 e-mail: info@wak-chemie.com • www.wak-chemie.com

### WAK - Chemie Medical GmbH



#### **Material Safety Data Sheet**

#### 1. **Product Description** Catalog No.: WAK-AQA-250-50L Product name: Acryl AquaClean Purpose of use: Microbiocidal additive for water soluble heating bath fluids for containers made of acrylic glass or aluminium. 2. **Chemical ingredients** Ketones, surfactants. 3. Health hazard data Weakly irritant for skin. Irritant for eyes or mucosa. Skin contact: promptly wash skin with water and mild soap. Eye or mucosa contact: rinse thoroughly with water for 15 minutes. Contact a physician. Wear appropriate protective clothing to avoid body contact. 4. **Emergency and first aid procedures** In case of ingestion: rinse mouth with water. Contact a physician. Never give anything by mouth to an unconscious person. Never make an unconscious person vomit. In case of allergic reactions contact a physician. 5. Fire and explosion hazard data Non inflammable. Water-soluble. 6. Steps to be taken in case material is released or spilled Use absorbent material for removal. Inactivate with 12% sodiumbisulfite-solution and rinse. 7. **Handling instructions** Avoid contact with eyes or mucosa. Wear protective clothing, gloves and glasses. 8. Storage Store light protected at room temperature. 9. **Physical data** 100°C **Boiling point:** Appearance: blue liquid pH-value: 6.8 Odor: n/a Specific weight: $1.01 \, a/cm^3$ water soluble Solubility: Flammable limits: none Chemical stability: stable 10. **Stability data** During decomposition by heat very small quantities of HCl, NOx und SO<sub>2</sub> may be set free. 11. **Toxicological data** Potentially toxic if ingested. 12. Waste disposal methods Diluted 1:200 liquid can be disposed like grey water. 13. **Transportation requirements** No special instructions. 14. General

While the information and recommendations set forth herein are believed to be accurate as of the date hereof, WAK-Chemie Medical GmbH makes no warranty with respect thereto and disclaims all liability from reliance thereon. Last revised 28. March 1999. 03/99-ConTaFree-1018ACA-SDB-E

Siemensstraße 9 D-61449 Steinbach / Ts. Germany Tel. + 49 61 71 / 2 84 30 • Fax + 49 61 71 / 69 56 46 e-mail: info@wak- chemie.com • www.wak- chemie.com

### WAK - Chemie Medical GmbH

